REDAX®
IT

KIT PERCUTANEO UNICO

versioni : "STANDARD e MULTI"

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INSTRUCTIONS FOR USE

"CENTESIS" GENERAL DESCRIPTION

The procedure involves the percutaneous insertion of a drainage catheter in order to remove fluids (liquids and air). All the versions are mainly but not exclusively dedicated to use in the patient's pleural cavity, abdomen, pericardium, etc. The purpose of this therapy is to alleviate the symptoms related to pneumothorax, even hypertensive, from different types of pleural effusions and deposits; the procedure may also involve the removal of liquids for diagnostic purposes or the infusion of therapeutic solutions.

INTENDED USE OF THE DEVICE

The "UNICO" system is mainly intended for creating a minimally invasive access to the various cavities of the patient so that accumulated liquids and air can be removed. The device can be used for gravity drainage or suction drainage, when connected to chest drainage systems, and also for obtaining fluid samples for diagnostic purposes. The device can also be used for the infusion of liquids such as physiological saline solutions for pleural lavage.

DESCRIPTION OF THE DEVICE

The "UNICO" device comes as a percutaneous kit that is available in different versions and various diameters. The kit is mainly composed of the following components:

1. A main body with three openings (1). The first opening (2), located at the rear, is for inserting the Verres needle or the atraumatic stylet. The two other openings (3) have Luer-lock connections and can be used for connecting drainage devices such as syringes, drainage bags, vacuum systems, etc. following the methods described on this instruction sheet. There is a selector (4) on the main body which, by rotating 90°, can be used to carry out either gravity drainage or the injection of liquids such as lavage solutions. The "basic" version has a drainage catheter (5) fixed to the main body and is equipped with holes and depth markings. In the version "with extension", the catheter is connected to the main body by means of a butterfly connector (6) which allows the positioning of an extension hose (7).
2. A Verres type needle (8) for introducing the catheter. This needle has a stylet with an atraumatic tip that protrudes from the end of the needle. This stylet can pass longitudinally until exposing the cutting edge of the needle during the insertion phase, which then returns to its extended position by means of a spring. The needle has a handle (9) with a visual indicator that informs the user of the position of the protective atraumatic stylet.
3. The devices are available in the "standard" version where the catheter is straight and in the "multi-use" version where the catheter has a pig-tail end.

ACCESSORIES

The kits normally contain a 2000 ml collection bag fitted with an air evacuation valve and an emptying valve (11). On request, other accessories can be provided:

- 60 ml syringe
- Single-use scalpel
- Drainage system connector
- Single-use clamp, etc.

The extension line (7) is included in the respective product versions.

WARNINGS AND GENERAL PRECAUTIONS

- **Read the Instructions for Use carefully before using the device;**
- **The use of this device is reserved for physicians and healthcare operators qualified to perform percutaneous access procedures and are aware of the risks and possible implications of these procedures.**
- **This device and all of its parts must only be used under safe conditions in the field of application and following the indications given on this instruction sheet. The manufacturer cannot be held responsible for any event deriving from inappropriate use or use that is different from that indicated.**
- **Check that the position indicator for the protective stylet moves freely and that the corresponding**


Fig. 8

indicator is green when the stylet is extended. If the green indicator does not appear during the initial compression on the outer surface of the thoracic wall, stop the insertion of the needle and replace the device.

- The device was designed for use with Redax products and accessories. If non Redax products and accessories are used, it is the physician's responsibility to check their compatibility before use.
- The device can be kept in situ for a maximum of 29 days. The length of the drainage procedure must be determined by the physician based on clinical considerations.
- Do not use the device if the outer packaging or any of its parts are damaged.
- Patients who have been diagnosed with a coagulopathy or other blood related illnesses must be carefully evaluated for the risk of bleeding.
- Do not manipulate the drainage in any way or make additional holes with cutting instruments.
- If fixing with a suture point, be extremely careful not to suture the hose or to puncture it with needles or cutting instruments. Avoid making the suture overly tight, as this could cause the drainage to stop.
- After the drain is positioned, ensure that the section with the drainage holes is completely inside the wound and that the connections are sealed perfectly.
- The drain must be removed manually by pulling with moderate pressure and avoiding any sudden movements. Avoid using metal instruments such as clamps or devices that could cause the drain to break.
- The drain may become difficult to remove if left in situ for long periods of time. More attention should be given during the removal operations.
- Once removed, ensure that the drain is completely intact. This is essential in order to rule out that any pieces of the drainage system have been left inside the cavity due to accidental lacerations.
- Single-use product. The reuse of this product could bring about changes in its performance and the risk of cross contamination.

PREPARATION FOR USE

Before using the system, carry out the preparation procedure as described below:

1. Check that the package is intact in order to ensure product function and sterility.
2. Using sterile procedure, open the packs making up the kit.
3. Determine the needle insertion position.
4. Create a sterile field around the preselected area and take the normal precautions for disinfecting the skin.
5. Remove the main part of the device with the connected catheter from the package; the catheter also comes pre-assembled to the body in the version with the butterfly connectors.
6. In the version with pig-tail, slide the rigid tube (Fig. 1) present on the catheter towards the end of the catheter until it is completely straightened out.
7. Insert the needle through the rear opening (2) until the handle is completely resting against the device. The needle and its protective stylet must protrude a few millimetres from the end of the catheter.
8. In the version with pig-tail, extract the rigid tube from the catheter.
9. Give a local anaesthetic at the site selected for the insertion of the needle using the normal anaesthetics.
10. If deemed necessary by the physician in order to help the insertion of the catheter, make a small incision in the skin.

WARNINGS

In cases of thoracentesis, the access point for needle insertion should be near the upper edge of the lower rib in order to avoid the neurovascular bundles located at the lower edge of each rib. For all other accesses, carefully check the position by means of radiographic

or ultrasound examination before insertion.

NEEDLE INSERTION AND CATHETER POSITIONING

After the previously described preparation has been performed, the needle can then be inserted.

WARNING: Even though the system is fitted with numerous safety devices, this operation is not without risks and it is therefore recommended that the procedure described below be followed.

1. Grip the device firmly with both hands. Use one hand on the main body to push in and guide the needle and hold its grip in contact with the device body. The other hand should be placed near the tip of the needle; gripping the catheter with the thumb and index finger, you can avoid inserting the needle too far into the patient's pleural cavity once you have gone through the skin.
2. Rest the tip of the needle against the patient's skin and apply light pressure until the green colour on the visual indicator disappears (Fig. 1).
3. Continue the insertion very slowly until the green on the indicator reappears; this means that the tip of the needle has entered the pleural space and the protective stylet is once again in the extended position, allowing its safe function (Fig. 2). The entrance into the cavity and the resulting reappearance of the green indicator may be accompanied by a clicking noise due to the activation of the spring that extends the protective stylet.
4. Move the device's handle in order to orient the insertion direction of the catheter. To check whether insertion has occurred in the desired spot, or rather in the area where there is the accumulation of fluid to be drained, an exploratory sample can be taken. To do this, connect the syringe to one of the two lower connectors, removing the corresponding cap. Draw the piston back and visually observe (Fig. 3) liquid entering the syringe; if this is the case, the syringe can be removed.
5. Begin extracting the needle and, at the same time, position the catheter, as shown in Fig. 4. In the version with pig-tail, when you remove the needle, the catheter will resume its rounded profile thanks to the shape memory of the material.

SUCTION OF FLUID (Pull-Push)

Once the catheter is positioned, the aspiration of the liquid in the pleural cavity can begin with the use of a syringe and the drainage bag. To connect these accessories, open the corresponding cap located on the two Luer-lock connections on the main body. Connect the syringe to one of the two connectors and the drainage bag to the other.

NOTE: the choice of which connector to use for the syringe or the bag is up to the physician, who can decide based on the ease of use. Draw the syringe piston (Fig. 5a); the liquid will be aspirated from the pleural cavity and into the syringe. By applying pressure to the piston, the liquid will automatically be sent to the drainage bag (Fig. 5b). By repeating this sequence of pulling and pushing, the liquid will be evacuated and collected in the bag.

NOTE: no additional operation is necessary and the selector must remain in its initial position (the knob is turned towards the catheter's axis, as seen in Fig. 5). The bag can be replaced or the syringe can be removed or both can be done at any time, depending on the result of the drainage. When removing accessories, it is recommended that the protective caps be closed in order to prevent possible contamination. If the caps are left open, there is however no risk of backflow towards the patient.

INSERTION OF THE EXTENSION HOSE (for the relative version)

The version of the device equipped with a butterfly connection is set up so an extension tube can be inserted between the catheter and the main body. As such, the drainage system can be easily fixed and held in place for long periods of time with greater patient comfort. To insert the extension tube, clamp the catheter with the disposable klemmer forceps included in the package. Now disconnect the butterfly connector and connect the tube as shown in Fig. 6. The klemmer forceps can be released after this operation. Proceed to fastening the catheter and connecting the accessories for fluid collection.


DRAINING OF FLUIDS: DRAINAGE BAG

After the catheter has been positioned and fixed in place, gravity drainage of the fluids can then begin. In this case, the drainage bag included in the package must be connected to one of the available connectors so that the draining liquids can be collected by means of gravity.

Proceed as follows for fluid drainage:

1. Connect the drainage bag to one of the Luer-lock connectors and ensure that the other connector is capped.
2. Position the drainage bag conveniently below the level of the patient.
3. Turn the selector on the main body until it stops, which corresponds to an angle of 90° with respect to the initial position

Fig. 7


Fig. 5b


Fig. 6

(Fig. 7). This represents a direct connection between the catheter and the drainage bag without the interpositioning of any valves.

WARNING: only use the drainage bag included in the kit or Redax replacement parts since they are fitted with special one-way valves that prevent backflow of fluids towards the patient.

- If the drainage bag is replaced, i.e. in case it is full, turn the selector to the initial position then replace the bag or the drain. After replacement, turn the selector 90°.

DRAINING OF FLUIDS: DRAINAGE SYSTEM

After positioning and securing the catheter, you can proceed with fluid drainage by gravity or suction using a dedicated postoperative drainage system (e.g. silicone reservoir, thoracic drainage system, etc.).

In this case, the drainage system must be connected to one of the available connectors using the special adaptor included in the package (12).

To connect a postoperative drainage device, operate as follows:

- Connect the adaptor included in the package to one of the Luer-lock connectors.
- Prepare the drainage system by carrying out the operations given on the manufacturer's instruction sheet for the respective products; then connect the system to the free end of the connector.
- Turn the selector on the main body until it stops, which corresponds to an angle of 90° with respect to the initial position. This represents a direct connection between the catheter and the drainage bag without the interpositioning of any valves.

WARNING: preferably use Redax products since they are fitted with special one-way valves that prevent the backflow of fluids towards the patient. If drainage devices from other manufacturer's are connected, make sure that a one-way, water seal or other type of valve is always present in order to ensure the correct evacuation of the pleural cavity.

- If the drainage system is to be replaced, i.e. in case it is full, turn the selector to the initial position then replace the bag or the drain. After replacement, turn the selector 90°.

DRAINING OF FLUIDS: SILICONE RESERVOIR FOR ASPIRATION

After the catheter has been positioned and fixed in place, the suction of fluids can begin using the specially designed silicone reservoir (optional).

Proceed as follows for fluid drainage:

- Connect the silicone reservoir to one of the Luer-lock connectors as shown in Fig. 8.
- Close both clamps located on the fluid inlet and on the bottom outlet towards the collection bag (if present). Now squeeze the pump until all of the air is evacuated from it. This operation is aided by the presence of a one-way valve.
- Close the protective cap on the one-way valve, as shown in Figure 8, and turn the selector on the main body until it stops, which corresponds to an angle of 90° with respect to the initial position. Open the clamp on the inlet of the reservoir and begin collecting the fluids.

WARNING: only use the Redax silicone reservoir, which is specially designed for use with the UNICO percutaneous drainage system.

- The silicone reservoir provides a constant vacuum of approximately 25cmH₂O until it is completely filled. When filled, the vacuum can be restored by simply connecting the special collection bag, opening the clamp and squeezing the reservoir until it is completely empty. Thanks to the valves included with the various devices, the vacuum is restored through emptying.

VOLUMETRIC MONITORING OF PATIENT'S AIR LEAKAGE

The volume of air leakage from the patient can be measured using the Redax silicone reservoir (optional) in the case of pleural cavity drainage. To perform this measurement, put the silicone reservoir under vacuum as described in the previous paragraph, evacuating all the air contained within. Open the clamp on the inlet of the reservoir and record the start time. The nursing staff can regularly check the degree of reservoir expansion until it is completely expanded. At this point, the volume of collected air is the reservoir's total capacity, unless there is any liquid collected, which must be subtracted. Below is a sample calculation:

Total volume of the silicone reservoir: 130 ml (including the dead space)
Time to total expansion: 20 minutes
Volume of collected liquid: 30 ml

Net volume of air collected: 130 – 30 = 100 ml
Patient air leakage = 100 ml / 20 min = 5 ml/min = 0.3 L/hr.

By evacuating the volume of fluid collected and restoring the reservoir's vacuum, the measurement can be repeated as many times as deemed necessary by the physician.

CONTRAINDICATIONS

Examination of medical literature indicates that, in very rare cases, the thoracentesis procedure can cause iatrogenic pneumothorax. This complication was related to manoeuvring errors when performing the thoracentesis and the use of high gauge needles and drainage systems.

There are no specific contraindications reported on the use of this device.

CLEANING AND MAINTENANCE

This is a single-use device and therefore requires no cleaning or maintenance.

DISPOSAL

After use, the device and all of its accessories may be contaminated with biological material and therefore must be disposed of as "hazardous biological" waste in accordance with the regulations in effect in the country where the device is discarded.

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